

University of Puerto Rico  
NMR LAB  
STANDARD OPERATING PROCEDURE: NMR-02

|                                 |                            |             |
|---------------------------------|----------------------------|-------------|
| TITLE: Proton $T_2$ Measurement |                            |             |
| Original Issue: 21 MAY 2006     | Revision Date: 21 AUG 2006 | Page 1 of 5 |

|                               |                               |
|-------------------------------|-------------------------------|
| Prepared By: José R. Martínez | Approved By: José R. Martínez |
|-------------------------------|-------------------------------|

### Introduction

Transverse relaxation ( $T_2$ ) is the mechanism by which the excited magnetization vector (conventionally shown in the  $x$ - $y$  plane) decays. This is always at least slightly faster than longitudinal relaxation. The magnitude of the magnetic moment in the  $x$ - $y$  plane decays according to

$$M = M_0 \exp\left(-\frac{t}{T_2}\right)$$


The CPMG (Carr-Purcell-Meiboom-Gill) ( $T_2$ ) experiment yields a signal of intensity

$$I_0 \exp\left(-\frac{t}{T_2}\right)$$


where  $t$  is the total evolution time ( $4n\Delta$ ). The value of  $\Delta$  ( $d_{20}$ ) in the pulse sequence should be much shorter than the reciprocal coupling constant  $1/J$  but long enough that the sample should not heat up significantly. A  $\Delta$  ( $d_{20}$ ) of 10 ms is usually appropriate. The experiment is repeated many times with different values of  $\tau$  and the resulting intensities used to find the value of  $T_2$ . This experiment works best for singlets. If the sample is very concentrated then the relaxation time will appear shorter than it really is due to saturation. In such a case, off-tune the probe, recalibrate the pulse widths and repeat the experiment.

We assume that the main program (TopSpin) to operate the instrument and the lock display ([lockdisp]-) were activated previously. It is also assumed that you previously ran a  $T_1$  experiment for the same sample.


**Figure 1:** Carr-Purcell-Meiboom-Gill Pulse Sequence














**Procedure**

| Step(s) and command(s)  | Comment(s)  |
|---|---|
| <p>1. [edc] ↵</p>  <p><b>Solvent:</b> the sample solvent<br/> <b>Experiment:</b> PROTON<br/> <b>Title:</b> T2 experiment of John Doe</p>   | <p>Create the data set to record a <sup>1</sup>H reference spectrum.</p> <p>↵ : means press enter in the keyboard.</p> <p><b>Important:</b><br/> DIR= /opt/topspin<br/> USER= NMR<br/> Experiment= PROTON</p>   |
| <p>In the command line (see Figure 2) type</p> <p>2. [ej] ↵ (to take out the previous sample)<br/> 3. introduce sample in spinner and measure sample depth<br/> 4. insert sample in magnet<br/> 5. [ij] ↵<br/> 6. [lock] ↵ → select sample solvent<br/> 7. [bsmsdisp] ↵ to adjust Z and Z2 or do [gradshim] ↵</p> | <p>Insert the sample in the magnet. Lock the spectrometer. Readjust the Z and Z2 shims until the lock level is optimized. If the instrument has gradients do gradient shimming. Tune and match the probehead for <sup>1</sup>H observation, if necessary.</p> |

*Figure 2: TopSpin Window*





| Step(s) and command(s)   | Comment(s)   |
|--|--|
|  <p>8. Type in the command line <b>[eda]</b> ↵, <b>PULPROG = zg</b>, type <b>solv</b> ↵ in the search window and press enter on keyboard; select the appropriate solvent and click on , or type in the command line <b>[getprosol]</b> ↵).</p>   | <p>You could type each parameter in the command line and set the value.</p> <p><b>[eda]</b> ≡ AcquPars</p> <p><b>PULPROG</b> ≡ Current pulse program</p> <p>The command  ≡ <b>[getprosol]</b> ↵, set the some acquisition parameters like pulse values.</p> |
| <p>9. Type in the command line <b>[d1] = 10</b> ↵, <b>[rga]</b> ↵ to adjust the receiver gain, then <b>[ns]=1</b> ↵, <b>[ds]=0</b> ↵ and <b>[zg]</b> ↵</p> <p>10. <b>[ef]</b>, <b>[apk]</b>; phase it if necessary.</p>  | <p><b>[d1]</b>: relaxation delay; 1-5*T<sub>1</sub></p> <p><b>[rga]</b>: Adjust receiver gain</p> <p><b>[ns]</b>: to select the number of acquisitions</p> <p><b>[ds]</b>: to dummy scans and</p> <p><b>[zg]</b>: to acquire the 1D proton spectrum</p>  |
| <p>11. Set <b>o1</b> on resonance with , then <b>[swh] = 500</b></p> <p>12. <b>[TD] = 1K</b> ; <b>[lb] = 2</b> ↵ ; <b>[ro]</b> ↵, then click on <b>Stop rotation</b> ; <b>[ns] = 8</b> ; <b>[ds] = 2</b> ; <b>[zg]</b></p> <p>13. <b>[ef]</b> ; <b>[apk]</b> ; confirm that no saturation, artifacts, etc. are observed.</p>  |    |
| <p>14. <b>[iexpno]</b> ↵</p> <p>15. <b>[eda]</b> ; <b>[pulprog] = cpmg</b> ; Switch from 1D to 2D clicking on  ; <b>[fnmode] = QF</b> ; <b>[aqmod] = qsim</b> or <b>DQD</b> (will depend on the instrument under use).</p> <p>16. <b>[td]</b> ↵ <b>F2 = 1024 F1 = 10</b></p> <p>17. <b>[ns] = 8</b> ; <b>[ds] = 2</b> ; <b>[d1] = 1-5 * T<sub>1</sub></b></p> <p>18. <b>[d20] = 10m</b> ; <b>[d11] = 30m</b> ; <b>[L4] = 10</b></p> | <p>From the previous data set create the new data set for the T<sub>2</sub> experiment. Also switch from 1D to 2D.</p> <p><b>[d1]</b>: relaxation delay; You should know the T<sub>1</sub> or an approximate value.</p> <p>The value for <b>L4</b> will depend on the values entered in the <b>vc</b> list</p>                                 |

|   |  |
|---|--|
| <p>19. [edlist] ↵, select <b>vc</b> (Variable counter) , select <b>t2delay</b> and click <b>Ok</b>. Enter the values as shown below:</p> <p>2<br/>20<br/>50<br/>100<br/>200<br/>300<br/>400<br/>500<br/>750<br/>1000</p> <p>Click <b>OK</b> (means save it)</p> <p>20. [zg]</p>   | <p><b>IMPORTANT:</b> The values shown are only an example; you have to set the values that better work for you.</p>  |
| <p>21. [edp] ↵<br/>         [si]↵ F2 = 512 F1 = 16<br/>         [wdw](F2) = EM<br/>         [lb] = 2<br/>         [phmod] = no<br/>         [pknl] = true<br/>         [bcmo] for F2 = quad, for F1 = no<br/>         [mc2] = QF</p>  | <p>[edp] ≡ ProcPars</p> <p>You could type each parameter in the command line and set the value.</p>  |
| <p>22. [xf2] ↵</p> <p>23. Click on  , and select at least two rows. Do a phase correction to that rows and store the correction (  ). Then .</p>   | <p>Save and return: </p> <p>Return: </p> |
| <p>24. Before any calculation open the home directory and in <b>Location</b> type the directory where is your <math>T_2</math> experiment file. For example,<br/> <b>Location:</b> file: /opt/topspin/data/user/nmr/NAME/EXPNO<br/>         In that Directory right click the mouse and select <b>Create New</b> → <b>Text File</b>, then name it <b>vdlist</b>. Open that file and enter the following list</p> <p>0.04<br/>0.4<br/>1<br/>2<br/>4<br/>6<br/>8<br/>10<br/>15<br/>20</p> | <p><b>user:</b> nmr in our case<br/> <b>NAME:</b> filemane of the experiment<br/> <b>EXPNO:</b> experiment number</p> <p><b>IMPORTANT:</b> The values shown are the delay values derived from the vclist.</p>      |

25. **Analysis** → **T1/T2 Relaxation** [t1guide]

26. Click on **Extract** → **Spectrum** → **Slice: 1**



27. Click on **Define Ranges** and select the signal () to which the  $T_2$  will be obtained.

28. Click on  to save as **Export Regions to Relaxation module and .ret**.

29. Click on **Relaxation Window**, select the **Fitting Function**  $\text{expdec: } I[t]=I[0] * \exp(-t/T)$  → **Function Type** → **expdec**. Confirm that the **vdlist** (List File Name) is selected too and click on **Apply** then **Ok**.

30. Select **Fitting Type** → **Area** and click on **Start Calculation**.

31. Click on **Display Report**, then **Print**.

#### References:

1. **TopSpin Users Guide**, Part Number H9469SA1 V2/April 1th 2004
2. 150\_and\_more\_v2.pdf.