



# **1D and 2D Experiments Step-by-Step Tutorial**

**Advanced Experiments  
User Guide**

**Version 003**



The information in this manual may be altered without notice.

BRUKER BIOSPIN accepts no responsibility for actions taken as a result of use of this manual. BRUKER BIOSPIN accepts no liability for any mistakes contained in the manual, leading to coincidental damage, whether during installation or operation of the instrument. Unauthorized reproduction of manual contents, without written permission from the publishers, or translation into another language, either in full or in part, is forbidden.

This manual was written by

Peter Ziegler

© February 15, 2008: Bruker Biospin Corporation

Billerica, Massachusetts, USA

P/N: B4722  
DWG-Nr.: 003

# Contents

|          |   |           |
|----------|---|-----------|
|          | <b>Contents .....</b>                                   | <b>3</b>  |
| <b>1</b> | <b>Introduction .....</b>                               | <b>7</b>  |
| 1.1      | General .....   | 7         |
| 1.2      | Disclaimer .....  | 7         |
| 1.3      | Warnings and Notes .....                                | 8         |
| 1.4      | Contact for Additional Technical Assistance .....       | 8         |
| <b>2</b> | <b>2-D Inverse Experiments .....</b>                    | <b>9</b>  |
| 2.1      | 2D edited HSQC .....                                    | 9         |
|          | Sample: .....   | 9         |
|          | Reference experiment .....                              | 9         |
|          | Setting up the HSQC experiment .....                    | 11        |
|          | Acquisition .....                                       | 11        |
|          | Processing .....  | 12        |
| 2.2      | 2D HMBC experiment .....                                | 13        |
|          | Sample: .....   | 13        |
|          | Preparation experiment .....                            | 13        |
|          | Setting up the HMBC experiment .....                    | 13        |
|          | Acquisition .....                                       | 14        |
|          | Processing .....  | 15        |
| 2.3      | Adding the F1 projection to the HSQC contour plot ..... | 16        |
|          | Creating the external projection spectrum .....         | 16        |
|          | Creating the plot template .....                        | 18        |
| <b>3</b> | <b>Diffusion Experiment .....</b>                       | <b>25</b> |
| 3.1      | Introduction .....                                      | 25        |
|          | Sample: .....   | 25        |
|          | Preparation experiment .....                            | 25        |
|          | Parameter set up .....                                  | 27        |
|          | Acquisition .....                                       | 29        |
|          | Processing .....  | 31        |
|          | Calculating the diffusion coefficient .....             | 33        |
|          | Displaying the DOSY plot .....                          | 39        |
| <b>4</b> | <b>Multiplet Analysis .....</b>                         | <b>41</b> |
| 4.1      | Multiplet assignments .....                             | 41        |
|          | Sample: .....   | 41        |
|          | Preparation experiment .....                            | 41        |
|          | Acquisition .....                                       | 42        |
|          | Processing .....  | 42        |
|          | Multiplet assignments .....                             | 45        |

|          |   |           |
|----------|---|-----------|
| <b>5</b> | <b>19F Experiments .....</b>  | <b>55</b> |
| 5.1      | Hardware necessary to observe 19F .....                                       | 55        |
|          | Probes .....  | 55        |
|          | Additional hardware .....   | 55        |
|          | 300 and 400MHz systems .....  | 55        |
|          | 500MHz and above .....  | 56        |
|          | Older AV systems .....  | 56        |
| 5.2      | 1-D 19F observe, without 1H decoupling .....                                  | 57        |
|          | Sample: .....   | 57        |
|          | Exploratory spectrum .....  | 57        |
|          | Optimizing the sweep width .....  | 59        |
|          | Baseline correction .....   | 59        |
| 5.3      | 19F observe, with 1H decoupling .....   | 61        |
|          | 1H observe, without 19F decoupling .....                                      | 62        |
|          | 1H observe, with 19F decoupling using WALTZ .....                             | 64        |
|          | 1H observe, with 19F decoupling using Garp .....                              | 66        |
|          | Plotting the 19F coupled and decoupled 1H spectrum on to<br>the same page. 66 |           |
| 5.4      | 2-D 19F/1H Hetero correlated experiment .....                                 | 68        |
|          | Sample: .....   | 68        |
|          | 1-D 19F reference spectrum .....  | 68        |
|          | 1-D 1H reference spectrum .....   | 70        |
|          | Setting up of the 2-D HETCOR experiment .....                                 | 71        |
| <b>6</b> | <b>1-D Selective NOESY .....</b>  | <b>75</b> |
| 6.1      | Introduction .....  | 75        |
|          | Reference spectrum .....  | 75        |
|          | Selective excitation region set up (example 1) .....                          | 77        |
|          | On resonance .....  | 77        |
|          | Setting up the acquisition parameters .....                                   | 78        |
|          | Running the experiment .....  | 78        |
|          | Selective excitation region set up (example 2) .....                          | 79        |
|          | Off resonance .....   | 79        |
|          | Setting up the acquisition parameters .....                                   | 81        |
|          | Running the experiment .....  | 81        |
|          | Selective excitation region set up (example 3) .....                          | 82        |
|          | Integration region file .....   | 82        |
|          | Calculating the selective pulse width and power level ....                    | 83        |
|          | Setting up the acquisition parameters .....                                   | 88        |
|          | Running the experiment .....  | 88        |
|          | Plotting the reference and the selective NOESY spectra on<br>the same page 89 |           |
| <b>7</b> | <b>1-D selective TOCSY .....</b>  | <b>91</b> |
| 7.1      | Introduction .....  | 91        |
|          | Reference spectrum .....  | 91        |
|          | Selective excitation region set up .....                                      | 93        |
|          | Off resonance .....   | 93        |
|          | Calculating the selective pulse width and power level ....                    | 94        |
|          | Setting up the acquisition parameters .....                                   | 98        |

|           |   |                   |
|-----------|---|-------------------|
|           | Running the experiment .....  | 98                |
|           | Plotting the reference and the TOCSY spectrum on to the same page.      | 99                |
| <b>8</b>  | <b><i>1-D DEPT using a shaped <sup>13</sup>C pulse .....</i></b>        | <b><i>101</i></b> |
| 8.1       | Introduction .....  | 101               |
|           | Experiment set up .....   | 101               |
|           | Calculating the shaped pulse power level .....                          | 102               |
|           | Running the experiment .....  | 105               |
| <b>9</b>  | <b><i>2-D HSQC using a shaped <sup>13</sup>C pulse .....</i></b>        | <b><i>107</i></b> |
| 9.1       | Introduction .....  | 107               |
|           | Reference spectrum .....  | 107               |
|           | Limit setting .....   | 108               |
|           | Running the 2-D HSQC using a 180 adiabatic inversion shaped pulse in F1 | 109               |
|           | Acquisition .....   | 110               |
|           | Processing .....  | 111               |
| <b>10</b> | <b><i>2-D Selective HMBC .....</i></b>                                  | <b><i>113</i></b> |
| 10.1      | Introduction .....  | 113               |
|           | Reference spectrum .....  | 113               |
|           | Limit setting .....   | 114               |
|           | Running a 2-D HMBC experiment .....                                     | 115               |
|           | Acquisition .....   | 116               |
|           | Processing .....  | 117               |
|           | Optimizing the parameters on the carbonyl region .....                  | 118               |
|           | Set up the selective pulse .....  | 119               |
|           | Setting up the acquisition parameters .....                             | 123               |
|           | Running the experiment .....  | 124               |



# *Introduction*

# 1

## *General*

1.1

---

This manual was written for AVANCE systems running TopSpin and should be used as a guide through the set up process for some experiments. The successful completion of the experiments in this manual presumes that all parameters have been entered in to the prosol table.

## *Disclaimer*

1.2

---

This guide should only be used for its intended purpose as described in this manual. Use of the manual for any purpose other than that for which it is intended is taken only at the users own risk and invalidates any and all manufacturer warranties.

Some parameter values, specially power levels suggested in this manual may not be suitable for all systems (e.g. Cryo probes) and could cause damage to the unit. Therefore only persons schooled in the operation of the AVANCE systems should operate the unit.

## Warnings and Notes

1.3

There are two types of information notices used in this manual. These notices highlight important information or warn the user of a potentially dangerous situation. The following notices will have the same level of importance throughout this manual.



---

Note: Indicates important information or helpful hints

---



---

**WARNING: Indicates the possibility of severe personal injury, loss of life or equipment damage if the instructions are not followed.**

---

## Contact for Additional Technical Assistance

1.4

For further technical assistance on the BPSU36-2 unit, please do not hesitate to contact your nearest BRUKER dealer or contact us directly at:

BRUKER BioSpin Corporation  
19 Fortune Drive, Manning Park  
Billerica, MA 01821  
USA

Phone: (978) 667-9580  
FAX: (978) 667-2955  
Email: [applab@bruker-biospin.com](mailto:applab@bruker-biospin.com)  
Internet: [www.bruker-biospin.com](http://www.bruker-biospin.com)



# 2-D Inverse Experiments

# 2

## 2D edited HSQC

2.1

### Sample:


30 mg Brucine in CDCl<sub>3</sub>

## Reference experiment

2.1.1

1. Run a 1D Proton spectrum, following the instructions in the Step-by-Step Tutorial, Basic Experiments User Guide, 1-D Proton Experiment, 2.2

Figure 2.1.




2. Type **wrpa 2** on the command line
3. Type **re 2**

- Expand the spectrum to display all peaks, leaving ca. 0.5 ppm of baseline on either side of the spectrum



**NOTE:** You may exclude the solvent peak, if it falls outside of the region of interest.

Figure 2.2.






- Click on  to set the sweep width and the O1 frequency of the displayed region

Figure 2.3.



- Write down the value of SW, rounding off to the nearest 1/10th of a ppm (e.g. **8.8** ppm)
- Write down the value of O1, rounding off to the nearest Hz (e.g. **1355** Hz)
- Click on 
- Type **sr** and write down the exact value

**Setting up the HSQC experiment**

2.1.2

1. Type **rpar HSQCEDETGP all**
2. Turn the spinner off




---

NOTE: 2-D experiments should be run non spinning

---

3. Select the '**AcquPars**' tab by clicking on it
4. Make the following changes:
 


|         |   |   |
|---------|---|---|
| SW [F2] | = | value from step 6, Preparation experiment 2.1.1 (e.g. <b>8.8</b> )  |
| O1 [Hz] | = | value from step 7, Preparation experiment 2.1.1 (e.g. <b>1355</b> ) |
| SOLVENT | = | <b>CDCl3</b>  |




---

All Bruker 2D inverse parameter sets use <sup>13</sup>C in the F1 dimension. Sweep width and O1 are optimized to include all Carbon peaks of interest. For HSQC and HMQC experiments the SW is optimized to 164 ppm.

---

5. Click on  to read in the Prosol parameters
6. Select the '**ProcPar**' tab by clicking on it
7. Make the following changes:
 

|         |   |  |
|---------|---|--|
| SR [F2] | = | value from step 9, Preparation experiment 2.1.1 (e.g. <b>0</b> ) |
|---------|---|--|
- 8 Select the '**Title**' tab by clicking on it
9. Change the title to: **2-D edited HSQC experiment of Brucine**
10. Select the '**Spectrum**' tab by clicking on it

**Acquisition**

2.1.3


1. Type **rga**
2. Type **zg** to start the acquisition



The standard Bruker parameter sets are optimized to run under complete automation through the use of AU programs. The name of the AU program is entered in the acquisition (eda) and processing (edp) parameter lists, as AUNM. For executing the processing AU program the command xaup may be used.

1. Type **edc2**

Figure 2.4.




2. Enter the EXPNO and PROCNO of the 1D Proton spectrum into the first column (data set 2)

3. Click on

4. Type **xaup**

Figure 2.5.





The processing AU program includes the 2D Fourier transform, phase correction, baseline correction and plotting of the data. The HSQC experiment is phase sensitive and it shows positive (red) peaks representing the CH and CH3 correlation and negative peaks (blue) shows the CH2.

## 2D HMBC experiment

2.2

### Sample:


20 mg Brucine in CDCl<sub>3</sub>

## Preparation experiment

2.2.1

1. Follow the instructions in the previous HSQC experiment 2.1.1 Preparation experiment step 1 through 9

Figure 2.6.



## Setting up the HMBC experiment

2.2.2

1. Type **rpar HMBCGP all**

2. Turn the spinner off



---

NOTE: 2-D experiments should be run non spinning

---

3. Select the '**AcquPars**' tab by clicking on it

4. Make the following changes:

SW [F2] = value from step 6, Preparation experiment 2.1.1 (e.g. **8.8**)

O1 [Hz] = value from step 7, Preparation experiment 2.1.1 (e.g. **1355**)


SOLVENT = **CDCI3**



---

All Bruker 2D inverse parameter sets use <sup>13</sup>C in the F1 dimension and the sweep width and O1 are optimized to include all Carbon peaks of interest. For HMBC experiments the SW is optimized to 220 ppm.

---

5. Click on  to read in the Prosol parameters

6. Select the '**ProcPar**' tab by clicking on it

7. Make the following changes:

SR [F2] = value from step 9 (Preparation experiment 2.1.1)

8 Select the '**Title**' tab by clicking on it

9. Change the title to: **2-D HMBC experiment of Brucine**

10. Select the '**Spectrum**' tab by clicking on it

---

### Acquisition

2.2.3

1. Type **rga**


2. Type **zg** to start the acquisition



The standard Bruker parameter sets are optimized to run under complete automation through the use of AU programs. The name of the AU program is entered in the acquisition (eda) and processing (edp) parameter lists, as AUNM. For executing the processing AU program the command xaup may be used.

1. Type **edc2**

Figure 2.7.




2. Enter the EXPNO and PROCNO of the 1D Proton spectrum into the first column (data set 2)

3. Click on

4. Type **xaup**

Figure 2.8.





The processing AU program includes the 2D Fourier transform, baseline correction and plotting of the data. The HMBC experiment uses magnitude mode for processing and shows only positive peaks.

### Adding the F1 projection to the HSQC contour plot

2.3



All Bruker 2D inverse parameter sets use the nucleus  $^{13}\text{C}$  in the F1 dimension. The sweep width and O1 frequency are set to include all Carbon peaks of interest. In most cases it is not necessary to run a Carbon spectrum to optimize the parameters. In the default plot template the F1 projection is disabled and therefore is not plotted.


If one wishes to add the F1 projection to the plot, then an additional  $^{13}\text{C}$  spectrum has to be obtained and a new plot template has to be created. For HMQC, HSQC type of experiments a DEPT45 or DEPT135 is best suited and for HMBC experiments a normal proton decoupled carbon spectrum should be used.

### Creating the external projection spectrum

2.3.1

1. Run a DEPT135 experiment following the instructions in the Step-by-Step Tutorial, Basic Experiments User Guide, DEPT135 Experiment 2.4.
2. Open the HSQC experiment
3. Type **edc2**

Figure 2.9.







4. Enter the EXPNO and PROCNO of the 1D Proton spectrum into the first column (data set 2)
5. Enter the EXPNO and PROCNO of the DEPT135 spectrum into the second column (data set 3)
6. Click on 
7. Type **xaup**
8. Right click inside the spectrum window

Figure 2.10.



9. Select '**Spectra Properties**' by clicking on it

Figure 2.11.






10. Enable '**Show projections**'
11. Enable Visible projections '**F1 + F2**'
12. Click on 
13. Right click on the F1 projection

Figure 2.12.



14. Select 'At Center' by clicking on it

Figure 2.13.




NOTE: Discard the plot

## Creating the plot template

2.3.2


1. Type **viewxwinplot**

Figure 2.14.



2. Click inside the spectrum window to display the green handles
3. Click on
4. Select the 'Basic' tab

Figure 2.15.




4. Make the following changes:

Position X = 3  
Dimension X = 15

- Click on 
- Select the **'2D Projection'** tab

Figure 2.16.




- Enable **'Data set left'**
- Make the following change:

Size = 3


- Click on 

Figure 2.17.



- Select the DEPT data set path
- Click on 
- Click on 
- Click on 

Figure 2.18.



15. Click on **1D/2D-Edit**

Figure 2.19.







Figure 2.21.




---

NOTE: Make sure to be in the directory [TopSpin home]\plot\layouts

---

25. Change Filename to **2D\_inv\_2proj.xwp**

26. Click on 





# Diffusion Experiment

# 3

## Introduction

3.1

---



---

NOTE: To run this experiment the instrument has to be equipped with the hardware to run Gradient experiments. Pulse field gradient NMR spectroscopy can be used to measure translational diffusion of molecules. The example in this chapter uses a mixture of two sugars dissolved in D<sub>2</sub>O.

---

### **Sample:**

Mixture of Glucose and Raffinose each 20mg in D<sub>2</sub>O


## Preparation experiment

3.1.1

---

1. Run a 1D Proton spectrum, following the instructions in the Step-by-Step Tutorial, Basic Experiments User Guide, 1-D Proton Experiment, 2.2

Figure 3.1.





2. Type **wrpa 2** on the command line
3. Type **re 2**
4. Expand the spectrum from 6ppm to 2ppm
5. Click on  to set the sweep width and the O1 frequency of the displayed region

Figure 3.2.





6. Click on 
7. Type **td 16k**
8. Type **si 8k**
9. Type **zg**
10. Type **ef**
11. Type **apk**
12. Type **abs**

Figure 3.3.



### Parameter set up

3.1.2




1. Type **ixpno**
2. Select the '**AcquPars**' tab by clicking on it
3. Click on  to display the pulse program parameters
4. make the following changes:
  - PULPROG = **stebpgp1s1d**
  - D20[s] = **0.1**
  - GPZ6[%] = **2**
  - GPZ7[%] = **-17.13**
  - P30[us] = **1800**
5. Type **rga**
6. Type **zg**
7. Type **ef**
8. Type **apk**
9. Type **abs**

Figure 3.4.



10. Type **ixpno**

11. Select the '**AcquPars**' tab by clicking on it

12. Click on  to display the pulse program parameters

13. make the following changes:

GPZ6[%] = **95**


14. Type **zg**


15. Type **ef**

16. Type **apk**

17. Type **abs**


Figure 3.5.




18. Click on  to open the multiple display window

19. Drag the previous experiment into the multiple display window (in this example it is experiment # 3) or type **re 3**

Figure 3.6.



NOTE: The intensity difference of the two spectra should be a factor of ~95. If the difference is less than 50, change P30 and or D20 in both data sets.

20. Click on  to exit the multiple display

## Acquisition

3.1.3



1. Type **ixpno**
2. Select the '**AcquPars**' tab by clicking on it
3. Make the following changes:  
PULPROG = **stebpgp1s**
4. Click on  to change the acqu dimension

Figure 3.7.



5. select '2D'

6. Click on


7. Make the following changes:

TDF1 = 16

FnMODE = QF


8. Type dosy on the command line

Figure 3.8.



9. Click on


Figure 3.9.



10. Enter 2 for first gradient amplitude

11. Click on


Figure 3.10.



12. Enter 95 for final gradient amplitude

13. Click on

Figure 3.11.



14. Enter **16** for the number of points

15. Click on

Figure 3.12.



16. Enter **l** for ramp type

17. Click on

Figure 3.13.





Figure 3.14.




18. Click on  to start the acquisition

## Processing

## 3.1.4

1. Select the **'Fid'** tab by clicking on it

Figure 3.15.



---

NOTE: Step 1 is only used to illustrate the DOSY experiment as a decay function.

---

2. Select the **ProcPars** tab by clicking on it

3. Make the following changes:

SI [F1] = 16

PH\_mod [F1] = no

PH\_mod [F2] = pk

4. Type **xf2** on the command line


5. Type **abs2** on the command line

6. Type **setdiffparm** on the command line

7. Select the **Spectrum** tab by clicking on it



Figure 3.16.



### Calculating the diffusion coefficient


3.1.5



NOTE: As you follow the steps below, message windows with important instructions will pop up. Please read this instructions very carefully.


1. Click on '**Analysis**' in the main menu
2. Select '**T1/T2 Relaxation**'

Figure 3.17.




3. Click on  to extract slice

Figure 3.18.



4. Click on


Figure 3.19.



5. Enter **1** for the slice number

6. Click on

Figure 3.20.





7. Click on  to define ranges

Figure 3.21.






8. Click on  Manual Integration


Figure 3.22.



9. Click on 

10. Define the regions by clicking the left mouse button and the use of the cursor lines

Figure 3.23.




11. Click on 

Figure 3.24.



12. Select **'Export Region To Relaxation Module'** by clicking on it



13. In the Guide window, click on  'Relaxation Window'


Figure 3.25.



14. Enable **'Area'**

15. In the guide window, click on  'Fitting Function'

Figure 3.26. 1






16. Click on 

Figure 3.27.




17. In the **'Fitting Function'** section, select **'vargrad'** and **'difflist'**

18. Click on 

19. In the guide window, click on  **'Start Calculation'**

Figure 3.28.





20. Click on 

Figure 3.29.




21. In the data window, click on  'Calculate fitting parameters for all data points'




NOTE: All calculated values are displayed in the 'Brief Report' section of the data window.

Figure 3.30.



22. In the guide window, click on  'Display Report'


Figure 3.31.



```

1 SIMFIT RESULTS
2 =====
3
4 Dataset : C:/data/pz/nmr/dosy/5/pdata/1/ct1t2.txt
5
6 AREA fit : Diffusion : Variable Gradient :
7
8 I=I[0]*exp(-D*SQR(2*PI*gamma*G*LD)*(BD-LD/3)*1e4)
9
10 16 points for Integral 1, Integral Region from 5.381 to 5.282 ppm
11
12 Converged after 61 iterations!
13
14 Results      Comp. 1
15
16 I[0]          = 1.499e-002
17 Diff Con.    = 3.539e-010 m2/s
18 Gamma        = 4.258e+003 Hz/G
19 Little Delta = 3.600m
20 Big Delta    = 99.900m
21
22 RSS          = 2.709e-008
23 SD           = 4.115e-005
24
25 Point      Gradient      Expt      Calc      Difference
26
27 1  6.740e-001  1.487e-002  1.497e-002  1.074e-004
28 2  2.765e+000  1.462e-002  1.463e-002  9.194e-006
29 3  4.855e+000  1.397e-002  1.389e-002  -7.354e-005
30 4  6.945e+000  1.283e-002  1.283e-002  -8.798e-006
31 5  9.036e+000  1.156e-002  1.151e-002  -4.687e-005
32 6  1.113e+001  1.009e-002  1.004e-002  -5.382e-005
33 7  1.322e+001  8.490e-003  8.514e-003  2.406e-005
34 8  1.531e+001  6.997e-003  7.019e-003  2.164e-005
35 9  1.740e+001  5.650e-003  5.625e-003  -2.583e-005
36 10 1.949e+001  4.347e-003  4.381e-003  3.378e-005
37 11 2.158e+001  3.297e-003  3.318e-003  2.044e-005
38 12 2.367e+001  2.412e-003  2.442e-003  3.033e-005
39 13 2.576e+001  1.752e-003  1.747e-003  -4.476e-006
40 14 2.785e+001  1.206e-003  1.215e-003  9.599e-006
41 15 2.994e+001  8.408e-004  8.216e-004  -1.919e-005
42 16 3.203e+001  5.552e-004  5.401e-004  -1.516e-005

```

23. In the guide window, click on  'Print Report

### Displaying the DOSY plot

3.1.6





1. Click on 
2. Click on  to close the NMR Relaxation Guide
3. Type **eddosy**

Figure 3.32.



4. Make the following change:

Scale = **Linear**

5. Click on  to set the upper and lower limits



6. Click on 

Figure 3.33.





# Multiplet Analysis

# 4

## Multiplet assignments

4.1

### Sample:


100 mg 2, 3,-Dibromopropionic acid in CDCl<sub>3</sub>

## Preparation experiment

4.1.1

1. Run a 1D Proton spectrum, following the instructions in the Step-by-Step Tutorial, Basic Experiments User Guide, 1-D Proton Experiment, 2.2
2. Type **ixpno** on the command line
3. Expand the spectrum to display all peaks, leaving ca. 0.5 ppm of baseline on either side of the spectrum

Figure 4.1.






4. Click on  to set the sweep width and the O1 frequency of the displayed region

Figure 4.2.



5. Click on 
6. Select the '**ProcPar**' tab by clicking on it
7. Make the following changes:  
LB = 0
- Select the '**Title**' tab by clicking on it
8. Change the title to: **1D Proton spectrum of 2, 3-Dibromopropionic acid**
9. Select the '**Spectrum**' tab by clicking on it

## Acquisition

4.1.2

1. Type **rga**
2. Type **zg** to start the acquisition

## Processing

4.1.3

1. Type **ft**
2. Type **apk**



---

**NOTE:** It may be necessary do a additional manual phase correction for a perfect phased spectrum.

---

3. Type **abs**



---

**NOTE:** If an internal reference such as TMS is added to the sample, a manual calibration should be done to the spectrum to assume a correct chemical shift of the

peaks. This may not be important for the multiplicity analysis, but for any spin simulation programs you may be using.



4. Expand the TMS peak
5. Click on  'Spectrum Calibration'
6. Move the cursor line into the center of the TMS peak
7. Click the left mouse button

Figure 4.3.






8. Change the value of the cursor frequency in ppm = 0
9. Click 
10. Expand the spectrum from 3.6 ppm to 4.6 ppm
11. Click with the right mouse button inside the spectrum window
12. Select 'Save Display Region To'
13. Enable the option 'Parameters F1/2 [dp1]'

Figure 4.4.



14. Click on **'Analysis'** in the main menu bar
  15. Select **'Peak Picking [pp]'** by clicking on it
- Figure 4.5.





16. Enable **'Define regions/peaks manually, adjust MI, MAXI'**
17. Click on 

Figure 4.6.






18. Move the cursor line to the left of the multiplet at 4.5 ppm
19. Click and hold the left mouse button and drag the cursor across the spectrum to the right of the multiplet at 3.7 ppm to draw a box over all multiplets
20. Click on  'Modify existing peak picking range'
21. Adjust the bottom line of the box to be above the baseline (Minimum intensity) and the top line above the highest peak of all multiplets (Maximum intensity)

Figure 4.7.



22. Click on 


**Multiplet assignments**

**4.1.4**

1. Expand the multiplet at 4.5 ppm

2. Click on '**Analysis**' in the main menu bar


Figure 4.8.




3. Select '**Multiplet Analysis**'

4. Select '**Enter Multiplet Mode [mana]**' by clicking on it


Figure 4.9.



5. Click on  'Define Multiplets Manually'


6. Place the cursor line to the left of the first peak of the multiplet

Figure 4.10.



7. Move the cursor line slowly towards the first peak
8. The cursor line will stop when it gets in to the center of the peak
9. Click the left mouse button

Figure 4.11.




10. Move the cursor line slowly towards the second peak
11. The cursor line will stop when it gets in to the center of the peak

12. Click the left mouse button



NOTE: A small marker is placed above the top of the first peak

Figure 4.12.




13. Move the cursor line in to the center of the two marked peaks

14. Click the right mouse button

15. Select 'Define Multiplet' by clicking on it


Figure 4.13.






16. Repeat steps 6 through 15 starting with the third peak and ending with the fourth peak

Figure 4.14.




17. Click on  'Couple Existing Multiplets'

18. Move the cursor line in to the center of the first two peaks marked 1

19. Click the left mouse button


Figure 4.15.



20. Move the cursor line in to the center of the second two lines marked 2

21. Click the left mouse button


Figure 4.16.



NOTE: While executing steps 20 through 21, the color of the brackets over the peaks 1 and 2 turn from black to red.

22. Move the cursor into the center of the displayed multiplet


Figure 4.17.



23. Click the right mouse button


24. Select 'Define Multiplet' by clicking on it

Figure 4.18.




25. Click the right mouse button inside the spectrum window
26. Select '**Finish Current Mode**' by clicking on it
27. expand the multiplet at 3.9 ppm
28. Repeat steps 6 through 26 for this multiplet

Figure 4.19.




29. Expand the multiplet at 3,7 ppm
30. Repeat steps 6 through 26 for this multiplet

Figure 4.20.



31. Display all 3 multiplets

Figure 4.21.





32. Click on  'Show Multiplet Report'

Figure 4.22.





33. Click on 

Figure 4.23.





34. Click on 

Figure 4.24.







---

NOTE: The connections are now assigned and the report can be printed.

---

35. Click on 

36. Click on  'Return, save multiplets [sret]'

# 19F Experiments

# 5

## Hardware necessary to observe 19F

5.1



NOTE: Below is a list of hardware options to observe or decoupled Fluorine on various Bruker systems and probes.

## Probes

5.1.1

|           |   |
|-----------|---|
| -QNP      | 19F/31P/13C/1H                          |
| -TXO      | 13C/1H/19F                              |
| -BBFO     | BB/19F/1H (300 and 400MHz systems only) |
| -BBO      | BB/1H (1H coil may be tunable to 19F)   |
| -BBI1H/BB | (1H coil may be tunable to 19F)         |
| -DUAL     | 1H/19F                                  |



NOTE: The probes listed above will have a Fluorine background with the exception of the Dual probe which is made Fluorine free. The BBO and BBI probes can only observe 19F without 1H decoupling. On the other hand, observing 13C and decoupling 19F is possible.

## Additional hardware

5.1.2

### 300 and 400MHz systems

-Internal amplifier      BLA-2BB

- 19F pass filter for doing observe 13C and 19F decoupling experiments.
- Other filters are built in to the preamplifiers (HPPR/2)



---

NOTE: By default amplifier 1 is connected to the X-BB preamplifier and amplifier 2 is connected to the 1H preamplifier. Each amplifier delivers 150 Watts from 10Mhz to the 31P frequency and 60 Watts above 31P to the 1H frequency and this will include the 19F frequency.

Standard pulse programs such as zg. zgdc etc. can be used to observe 19F.

---

### **500MHz and above**

- external amplifiers      BLAXH (less then 1.5 years old)  
   BLA(R)H, BLAX combinations
- external QNP accessory unit for RF routing
- 19F pass filter for doing observe 13C and 19F decoupling experiments.
- Other filters are built in to the preamplifiers (HPPR/2)



---

NOTE: The 19F signal is generated on the 1H stage of the amplifier and the QNP accessory unit is designed to route the 19F frequency either to the 19F selective or X-QNP output. In addition it switches between the 1H and 19F for decoupling either of the nuclei.

Pulse programs have to include the routing and switching statements such as QNP\_X, QNP\_F, SWITO\_F, SWITO\_H.

---

### **Older AV systems**

- external amplifier      BLAXH (more then 1.5 years old)
- The QNP switch unit is built in to the amplifier and the functions are the same as the above QNP accessory unit.
- 19F pass filter for doing observe 13C and 19F decoupling experiments.



-Additional filters such as 'Band Pass X, 19F//Band Stop 1H' and 'Band Pass 1H//Band Stop 19F' are necessary if a HPPR/1 is in use.



NOTE: The 19F signal is generated on the 1H stage of the amplifier and the QNP accessory unit is designed to route the 19F frequency either to the 19F selective or X-QNP output. In addition it switches between the 1H and 19F for decoupling either of the nuclei.

Pulse programs have to include the routing and switching statements such as QNP\_X, QNP\_F, SWITO\_F, SWITO\_H.

## 1-D 19F observe, without 1H decoupling

5.2



The 19F chemical shift range is rather large and covers approximately from +100ppm to -300ppm. The default sweep width of the Bruker standard 19F parameter sets may not cover the whole chemical shift range and adjustment may be needed. A common reference standard is: CFCI<sub>3</sub> at 0ppm. Others standards such as CF<sub>3</sub>COOH and C<sub>6</sub>F<sub>6</sub> may also be used.

### Sample:

2,2,3,4,4,4-Hexafluoro-1-butanol in Acetone-d<sub>6</sub>

CF<sub>3</sub>-CFH-CF<sub>2</sub>-CH<sub>2</sub>-OH

## Exploratory spectrum

5.2.1

1. Type **edc** and change the following parameters

Figure 5.1.

Prepare for a new experiment by creating a new data set and initializing its NMR parameters according to the selected experiment type. For multi-receiver experiments several datasets are created. Please define the number of receivers in the box below.

NAME: f19exp  
 EXPNO: 1  
 PROCNO: 1  
 DIR: C:  
 USER: pz  
 Solvent: Acetone  
 Experiment Dirs.: C:/Bruker/TOPSPIN2.1.1/exp/star/nmr/par  
 Experiment: F19  
 TITLE: 1-D 19F observe, without 1H decoupling  
 2,2,3,3,4,4-Hexafluoro-1-butanol  
 1 Receivers (1,2, ...8)

OK Cancel More Info... Help




2. Click on 
3. Insert the sample
4. Type **lock** and select Acetone
5. Shim for best homogeneity
6. Tune the probe
7. Select the '**AcquPars**' tab by clicking on it
8. Click on  to read in the Prosol parameters
9. Type **rga**
10. Type **zg** to start the acquisition
11. Process and Phase correct the spectrum

Figure 5.2.






In this example, the right most peak at ca. 220ppm is too close to the edge and may be distorted by the digital filtering. In this case, the SW and O1P should be adjusted.


1. Select the '**AcquPars**' tab by clicking on it
2. Change the following parameters:  
 SW [PPM] = **200**  
 O1P [PPM] = **-140**
3. Type **zg** to start the acquisition
4. Process and Phase correct the spectrum

Figure 5.3.



1. Display the full spectrum
2. Expand the spectrum vertically

Figure 5.4.



If a Fluorine background signal is present, a simple abs will not straighten the baseline and a linear prediction calculation may be necessary. See steps below.

3. Type **convdta**

Figure 5.5.



4. Type **2** into the convdta window

5. Click on

6. Select the '**Procpar**' tab by clicking on it

7. Change the following parameters:

ME\_mod = **LPbc**

NCOEF = **32**


TDoff = **16**

8. Type **ef**

9. Phase correct the spectrum

10. Type **abs**

Figure 5.6.



## 19F observe, with 1H decoupling

5.3

1. Type **edc** and change the following parameters

Figure 5.7.

2. Click on **OK**

3. Insert the sample

4. Type **lock** and select Acetone

5. Shim for best homogeneity



6. Tune the probe
7. Select the '**AcquPars**' tab by clicking on it
8. Change the following parameters:
  - SW [PPM] = **200**
  - O1P [PPM] = **-140**
9. Click on  to read in the Prosol parameters
10. Select the '**Spectrum**' tab by clicking on it
11. Type **rga**
12. Type **zg** to start the acquisition
13. Process and Phase correct the spectrum
14. To get rid of the background signal, follow the instructions in 5.2.3, Baseline correction, steps 1 through 9

Figure 5.8.

**1H observe, without 19F decoupling****5.3.1**

1. Type **edc** and change the following parameters

Figure 5.9.

Prepare for a new experiment by creating a new data set and initializing its NMR parameters according to the selected experiment type. For multi-receiver experiments several datasets are created. Please define the number of receivers in the box below.

NAME: f19exp  
 EXPNO: 1  
 PROCNO: 1  
 DIR: C:  
 USER: pz  
 Solvent: Acetone  
 Experiment Dirs.: C:/Bruker/TOPSPIN2.1/exp/stan/nmr/par  
 Experiment: PROTON  
 TITLE: 1-D 1H observe, without 19F decoupling  
 2,2,3,4,4,4-Hexafluoro-1-butanol  
 1 Receivers (1,2, ...8)

OK Cancel More Info... Help




2. Click on 
3. Insert the sample
4. Type **lock** and select Acetone
5. Shim for best homogeneity
6. Tune the probe
7. Select the '**AcquPars**' tab by clicking on it
8. Click on  to read in the Prosol parameters
9. Select the '**Spectrum**' tab by clicking on it
10. Type **rga**
11. Type **zg** to start the acquisition
12. Process and Phase correct the spectrum

Figure 5.10.



### 1H observe, with 19F decoupling using WALTZ

5.3.2

1. Type **ixpno**




NOTE: The experiment number is incremented by one (e.g. EXPNO = 2)

2. Type **rpar PROF19DEC all**
3. Tune the probe for 19F and 1H
4. Select the '**AcquPars**' tab by clicking on it
5. Change the following parameters:

TD = **64k**  
 DS = **10**  
 O2P [PPM] = **-180**  
 SOLVENT = **Acetone**

6. Select the '**ProcPars**' tab by clicking on it
- Change the following parameter:


SI = **32k**

7. Click on  to read in the Prosol parameters
8. Select the '**Title**' tab by clicking on it



9. Change the title to: **1-D 1H experiment with 19F decoupling using WALTZ16 2,2,3,4,4,4-Hexafluoro-1-Butanol**
10. Select the 'Spectrum' tab by clicking on it
11. Type **rga**
12. Type **zg** to start the acquisition
13. Process and Phase correct the spectrum

Figure 5.11.




---


The Bruker standard parameter set PROF19DEC is using WALTZ for decoupling  $^{19}\text{F}$ . This may not be sufficient of a bandwidth to cover the  $^{19}\text{F}$  chemical shift range of some of the  $^{19}\text{F}$  spectra. In this example the  $^{19}\text{F}$  signals covers a sweep width of 200 ppm. To decouple all the  $^{19}\text{F}$  peaks, two approaches can be applied. Using the WALTZ decoupling the O2 frequency would have to be adjusted for the various  $^{19}\text{F}$  resonances which results in multiple proton spectra. Using garp or adiabatic pulses widens the decoupling range. Below is a example using garp decoupling.

---

1. Type **ixpno** in the command line




NOTE: The experiment number is incremented by one (e.g. EXPNO = 3)

2. Select the '**AcquPars**' tab by clicking on it
3. Click on  to display the pulse program parameters
4. Make the following changes:
 


|         |   |  |
|---------|---|--|
| CPDPRG2 | = | <b>garp</b>  |
| PCPD2   | = | <b>70</b>  |
| PI12    | = | <b>calculate the power level in the prosol table</b> |
5. Type **rga**
6. Type **zg** to start the acquisition
7. Process and Phase correct the spectrum

Figure 5.12.



Plotting the 19F coupled and decoupled 1H spectrum on to the same page.

5.3.4

1. Type **re 3** to display the 1H observe with 19F decoupled spectrum
- 2.. Click on 




3. Type **re 1** on the command line (1H observe without 19F decoupled spectrum)
4. Click on  to separate the two spectra
5. Using the display tools  to adjust the spectra

Figure 5.13.



6. Type **print** on the command line to print the active window




---

NOTE: To plot the two spectra using the plot editor, follow the instructions in the manual Step-by-Step Tutorial, Basic Experiments Users Guide, Chapter 7, Homo-decoupling, 7.1.6 Plotting the reference and decoupled spectra on the same page, steps 1 through 21.

---



There are currently no standard parameter sets for 19F 2-D experiments. The instructions below will guide you through the creation of the 19F 2-D parameter set and running the experiment.

**Sample:**

2,2,3,4,4,4-Hexafluoro-1-butanol in Acetone-d6  
 $\text{CF}_3\text{-CFH-CF}_2\text{-CH}_2\text{-OH}$

1. Type **edc** and change the following parameters

Figure 5.14.

2. Click on **OK**
3. Insert the sample
4. Type **lock** and select Acetone
5. Shim for best homogeneity
6. Tune the probe
7. Select the '**AcquPars**' tab by clicking on it




8. Click on  to read in the Prosol parameters
9. Select the 'Spectrum' tab by clicking on it
10. Type **rga**
11. Type **zg** to start the acquisition
12. Process and Phase correct the spectrum
13. Expand the spectrum to display all peaks, leaving ca. 15ppm of baseline on either side of the spectrum
14. Click on  to set the sweep width and the O1 frequency of the displayed region

Figure 5.15.





15. Write down the SW value, rounding off to the nearest ppm (e.g. **170**)
16. Write down the O1 value, rounding off to the nearest Hz (e.g. **-40944**)
17. Click on 
18. Type **sr** and write down the exact value (e.g. **0**)
19. Type **zg** to start the acquisition
20. Process and Phase correct the spectrum

Figure 5.16.



1. Type **ixpno**



NOTE: The experiment number is incremented by one (e.g. EXPNO = 2)




2. Type **rpar PROTON all**
3. Select the 'AcquPars' tab by clicking on it
4. Click on  to read in the Prosol parameters
5. Select the 'Spectrum' tab by clicking on it
6. Type **rga**
7. Type **zg** to start the acquisition
8. Process and Phase correct the spectrum
9. Expand the spectrum to display all peaks, leaving ca. 0.5 ppm of baseline on either side of the spectrum
10. Click on  to set the sweep width and the O1 frequency of the displayed region

Figure 5.17.





11. Write down the SW value, rounding off to the nearest 1/10th of a ppm (e.g. **3**)
12. Write down the O1 value, rounding off to the nearest Hz (e.g. **1320**)
13. Click on 
14. Type **sr** and write down the exact value (e.g. **0**)
15. Type **zg** to start the acquisition
16. Process and Phase correct the spectrum

Figure 5.18.



### Setting up of the 2-D HETCOR experiment

5.4.3

1. Type **iexpno**



NOTE: The experiment number is incremented by one (e.g. EXPNO = 3)

2. Type **rpar HCCOSW all**
3. Turn the spinner off




NOTE: 2-D experiments should be run non spinning

4. Type **edasp**
5. Make the following change:  

$$\text{NUC1} = 19\text{F}$$

Figure 5.19.



6. Click on 

7. Select the '**AcquPars**' tab by clicking on it

8. Make the following changes:

PULPROG = **hfcoqfqm**


SW F2 [ppm] = e.g. **170** (19F reference spectrum 5.4.1, step 15)


SW F1 [ppm] = e.g. **3** (1H reference spectrum 5.4.2, step 11)

O1 [Hz] = e.g. **-40944** (19F reference spectrum 5.4.1, step 16)

O2 [Hz] = e.g. **1320** (1H reference spectrum 5.4.2, step 12)

SOLVENT = **Acetone**

9. Click on  to read in the Prosol parameters

10. Click on  to display the pulse program parameters

11. Make the following changes:

CNST2 = **25** = J(FH)

12. Select the '**ProcPar**' tab by clicking on it

13. Make the following changes:

SR F2 = e.g. **0** (19F reference spectrum 5.4.1, step 18)

SR F1 = e.g. **0** (1H reference spectrum 5.4.2, step 14)

WDW F2 = **SINE**

WDW F1 = **SINE**

SSB F2 = **2**

SSB F1 = **2**

14. Select the '**Title**' tab by clicking on it

15. Change the title to: **2-D 1H/19F HETCOR experiment**  
**2,2,3,4,4,4-Hexafluoro-1-Butanol**

16. Select the '**Spectrum**' tab by clicking on it


17. Type **rga**

18. Type **zg** to start the acquisition



19. Type **xfb**
20. Select the reference spectra as the external projections
21. Adjust the contour level

Figure 5.20.





# 1-D Selective NOESY

# 6

## Introduction

6.1



NOTE: To run this experiment the instrument has to be equipped with the hardware to do Shaped Pulses and Gradients. Three different ways to run this experiment are discussed in this chapter and can also be applied to other selective experiments such as SELCOSY, SELROESY and SELTOCSY.

### Sample:

30 mg Paomic acid in DMSO

## Reference spectrum

6.1.1


1. Type **edc** and change the following parameters

Figure 6.1.

Prepare for a new experiment by creating a new data set and initializing its NMR parameters according to the selected experiment type. For multi-receiver experiments several datasets are created. Please define the number of receivers in the box below.

|                       |  |
|-----------------------|--|
| NAME                  | selexp                                     |
| EXPNO                 | 1  |
| PROCNO                | 1  |
| DIR                   | C:   |
| USER                  | pz   |
| Solvent               | DMSO                                       |
| Experiment Dirs.      | C:/Bruker/TOPSPIN2.1.1/exp/stan/nmr/par    |
| Experiment            | PROTON                                     |
| TITLE                 | 1-D Proton experiment<br>30 mg Pamoic acid |
| Receivers (1,2, ...8) | 1  |

OK Cancel More Info... Help

2. Click on 
3. Insert the sample
4. Type **lock** and select DMSO
5. Turn the spinner off



---

NOTE: selective excitation experiments should be run non spinning

---



6. Shim for best homogeneity
7. Select the '**AcquPars**' tab by clicking on it
8. Click on  to read in the Prosol parameters
9. Tune the probe
10. Type **rga**
11. Process and Phase correct the spectrum

Figure 6.2.



**On resonance**

NOTE: Make sure that the SW is large enough to cover the entire Spectrum accounting for the position of O1. The shaped pulse is applied on resonance (at the o1 position) The power level and width of the excitation pulse have to be known and entered into the Prosol parameter table




1. Type **wrpa 2**
2. Type **re 2**
3. Select the **'Title'** tab by clicking on it
4. Change the title to: **Selective NOESY experiment**
5. Select the **'Spectrum'** tab by clicking on it
6. Expand the signal region at 8.5 ppm
7. Click on 


Figure 6.3.



8. Move the cursor line to the center of the peak and click the left mouse button


Figure 6.4.



9. Click on 

### Setting up the acquisition parameters

6.1.3

1. Select the '**AcquPars**' tab by clicking on it
2. Click on  to display the pulse program parameters
3. Make the following changes:

PULPROG = **selnogg**  
 NS = **64**  
 DS = **8**  
 D1 = **2**  
 D8 = **0.700**  
 SPNAM2 = **Gaus1.1000**  
 SPOFF2 = **0**  
 GPNAM1 = **sine.100**  
 GPNAM2 = **sine.100**  
 GPZ1 = **15**  
 GPZ2 = **40**



**NOTE:** The mixing time D8 is dependent on the size of the Molecule and the magnetic strength. It can vary from a large Molecule to a small one from 100 ms to 800 ms.

### Running the experiment

6.1.4


1. Select the '**Spectrum**' tab by clicking on it
3. Type **zg** to start the acquisition
4. Type **ef**

## 5. Phase the spectrum using the manual phase adjust



NOTE: Phase the selective excited peak negative to be sure the correct phase of the noe peaks.

Figure 6.5.

**Selective excitation region set up (example 2)****6.1.5****Off resonance**

NOTE: This method does not require a large SW. The shaped pulse is applied off resonance (not on the O1 position). The power level and pulse width of the excitation pulse have to be known and entered into the Prosol parameters.

1. Run a Reference spectrum, following the instructions in 6.1.1 Reference Spectrum in this Chapter.
2. Type **wrpa 2**




3. Type **re 2**
4. Select the **'Title'** tab by clicking on it
5. Change the title to: **Selective NOESY experiment**
6. Select the **'Spectrum'** tab by clicking on it
7. Expand the signal region at 8.5 ppm
8. Click on 


Figure 6.6.




9. Move the cursor line to the center of the peak and click the left mouse button

Figure 6.7.



10. Write down the O1/2/3 (Hz) value showing in the Info window (e.g. **2542.9**)
12. Click on 
13. Type **O1** and write down the current value (e.g. **1853.43**)
14. Calculate the difference of step 9 and 11 (e.g. **689.47**)




15. Click on 



NOTE: If the signal is down field of O1, a positive value must be entered for spoff. If the signal is up field of O1, spoff will have a negative value.

### Setting up the acquisition parameters

6.1.6

1. Select the '**AcquPars**' tab by clicking on it
2. Click on  to display the pulse program parameters
3. Make the following changes:

PULPROG = **selnogp**  
 NS = **64**  
 DS = **8**  
 D1 = **2**  
 D8 = **0.700**  
 SPNAM2 = **Gaus1.1000**  
 SPOFF2 = value from 6.1.5, step 14 (e.g. **689.47**)  
 GPNAM1 = **sine.100**  
 GPNAM2 = **sine.100**  
 GPZ1 = **15**  
 GPZ2 = **40**



NOTE: The mixing time D8 is dependent on the size of the Molecule and the magnetic strength. It can vary from a large Molecule to a small one from 100 ms to 800 ms.

### Running the experiment

6.1.7


1. Select the '**Spectrum**' tab by clicking on it

2. Type **zg** to start the acquisition
3. Type **ef**
4. Phase the spectrum using the manual phase adjust



NOTE: Phase the selective excited peak negative to be sure the correct phase of the noe peaks.

Figure 6.8.



## Selective excitation region set up (example 3)

6.1.8

### Integration region file



NOTE: In this example the shaped pulse is applied at a position determined using an integration region file and therefore does not require a large SW. This method calculates the precise shaped pulse for the selected peak using the 90 degree hard pulse and the Shape Tool program.

1. Run a Reference spectrum, following the instructions in 6.1.1 Reference Spectrum in this Chapter.
2. Type **wrpa 2**
3. Type **re 2**
4. Select the 'Title' tab by clicking on it
5. Change the title to: **Selective NOESY experiment**
6. Select the 'Spectrum' tab by clicking on it
7. Type **pulprog selnogg** in the command line

### Calculating the selective pulse width and power level


6.1.9



In this example the shaped pulse width and power level are determine using the '**Calculate the excitation profiles**' option in the shaped tool program. Other method of calculating the pulse width and power level can be used, see Chapter 7, 1-D Selective TOCSY, Bandwidth region file, in this manual, or use the Prosol parameters to run this experiment.


1. In the main menu click on '**Spectrometer**' and select '**Shape Tool**' or type **stdisp** in the command line
2. In the shape tool menu bar click on  and select '**Shape**'

Figure 6.9.




3. Select '**Gaus1.1000**'
4. Click on 

Figure 6.10.



5. In the shapetool window expand the signal region around 8.5 ppm


Figure 6.11.



6. Click on  'New excitation region'

7. Adjust the excitation region to cover the full width of the peak

Figure 6.12.





8. Click on  'Calculate the excitation profiles'

Figure 6.13.



9. In the main menu click on 'Options' and select 'Define Parameter Table' by clicking on it

Figure 6.14.





Figure 6.15.




10. Make the following changes:

- Length of shaped pulse = **p12**
- Power Level of shaped pulse = **SP2**
- Name of shaped pulse = **SPNAM2**


11. Click on

12. Click on

Figure 6.16.




13. Select a new name or use the default name displayed

14. Click on 



NOTE: If the shape region file already exists, an overwrite message will appear like the one below.

Figure 6.17.







15. Click on 


Figure 6.18.



NOTE: The Target data set window above is to verify the correct data set and can be switched off by enable the 'Do not ask again' option.

16. Click on 

17. Click on  to close the Shape Tool window

1. Select the '**AcquPars**' tab by clicking on it
2. Click on  to display the pulse program parameters
3. Make the following changes:

NS = 64  
DS = 8  
D1 = 2  
D8 = 0.700  
GPNAM1 = sine.100  
GPNAM2 = sine.100  
GPZ1 = 15  
GPZ2 = 40



---

**NOTE:** The mixing time D8 is dependent on the size of the Molecule and the magnetic strength. It can vary from a large Molecule to a small one from 100 ms to 800 ms.

---

1. Select the '**Spectrum**' tab by clicking on it
2. Type **zg** to start the acquisition
3. Type **ef**
4. Phase the spectrum using the manual phase adjust




---

**NOTE:** Phase the selective excited peak negative to a sure the correct phase of the noe peaks.

---



Figure 6.19.




### Plotting the reference and the selective NOESY spectra on the same page

6.1.12

1. Type **re 2** to display the selective NOESY spectrum
2. Click on
3. Type **re 1** on the command line (reference spectrum)
4. Click on
5. Using the display tools to adjust the spectra

Figure 6.20.



6. Type **prnt** on the command line to print the active window



---

NOTE: To plot the two spectra using the plot editor, follow the instructions in the manual Step-by-Step Tutorial, Basic Experiments Users Guide, Chapter 7, Homodecoupling, 7.1.6 Plotting the reference and decoupled spectra on the same page, steps 1 through 21.

---

# 1-D selective TOCSY

# 7

## Introduction

7.1



NOTE: To run this experiment the instrument has to be equipped with the hardware to do Shaped Pulses and Gradients. The method to determine the pulse width and power level for the selective pulse in this chapter, can also be used for other selective experiments such as SELCOSY, SELROESY and SELNOESY.

### Sample:

50 mM Gramicidin S in DMSO

## Reference spectrum

7.1.1


1. Type **edc** and change the following parameters

Figure 7.1.

Prepare for a new experiment by creating a new data set and initializing its NMR parameters according to the selected experiment type. For multi-receiver experiments several datasets are created. Please define the number of receivers in the box below.

|                       |   |
|-----------------------|---|
| NAME                  | seitocsy  |
| EXPNO                 | 1   |
| PROCNO                | 1   |
| DIR                   | C:  |
| USER                  | pz  |
| Solvent               | DMSO  |
| Experiment Dirs.      | C:/Bruker/TOPSPIN2.1.1/exp/stan/nmr/par             |
| Experiment            | PROTON  |
| TITLE                 | 1-D Proton spectrum<br>50mM Gramicidin S in DMSO-d6 |
| Receivers (1,2, ...8) | 1   |

OK Cancel More Info... Help

2. Click on 
3. Insert the sample
4. Type **lock** and select DMSO
5. Turn the spinner off



NOTE: selective excitation experiments should be run non spinning



6. Shim for best homogeneity
7. Select the '**AcquPars**' tab by clicking on it
8. Click on  to read in the Prosol parameters
9. Tune the probe
10. Type **rga**
11. Process and Phase correct the spectrum

Figure 7.2.



**Off resonance**

NOTE: In this example the shaped pulse is applied at the off resonance position and therefore does not require a large SW. Other excitation region set up method can be used to run this experiment, see Chapter 2, 1-D Selective NOESY in this manual.




1. Type **wrpa 2**
2. Type **re 2**
3. Select the **'Title'** tab by clicking on it
4. Change the title to: **Selective TOCSY experiment**
5. Select the **'Spectrum'** tab by clicking on it
6. Expand the amid peak of Leucine at 8.3 ppm
7. Click on 



Figure 7.3.



8. Move the cursor line to the center of the peak and click the left mouse button

Figure 7.4.



9. Write down the O1/2/3 (Hz) value showing in the Info window (e.g. 2498.9)
10. Click on 
11. Type **O1** and write down the current value (e.g. 1853.43)
12. Calculate the difference of step 9 and 11 and write down the value, (e.g. 645.47 Hz)
14. Click on 




---

NOTE: If the signal is down field of O1, a positive value must be entered for spoff. If the signal is up field of O1, spoff will have a negative value.

---

### Calculating the selective pulse width and power level

7.1.3




---

In this example the shaped pulse width and power level are determine using the '**Calculate Bandwidth**' option in the shaped tool program. Other method of calculating the pulse width and power level can be used, see Chapter 2, 1-D Selective NOESY, integration region file, in this manual, or use the Prosol parameters to run this experiment.

---



1. Click on 
2. Position the cursor line at the left side of the peak, up 1/5 from the baseline
3. Click the left mouse button and drag the cursor line to the right side of the peak, up 1/5 from the baseline

Figure 7.5.





4. Write down the value for the distance between the two cursor lines (e.g. 19)
5. Type **pulprog selmlgp**
6. In the main menu click on **'Spectrometer'** and select **'Shape Tool'** or type **stdisp**
7. In the shape tool menu bar click on  and select **'Open Shape'**

Figure 7.6.






8. Select **'Gaus1\_180r.1000'**
9. Click on 
10. In the main menu click on **'Analysis'** and select **'Calculate Bandwidth for Refocusing -My'**

Figure 7.7.



11. Type the value from step 4 (e.g. 19) in to the Calculator window 'Delta Omega [Hz]' and hit the Enter key


Figure 7.8.



NOTE: The value for 'Delta T [usec]' is calculated after executing step 11.

12. In the main menu click on 'Options' and select 'Define Parameter Table'

Figure 7.9.




13. Make the following changes:


Length of shaped pulse = p12



Power Level of shaped pulse = **SP2**  
 Name of shaped pulse = **SPNAM2**

Figure 7.10.



14. Click on 





15. Click on 


Figure 7.11.



**NOTE:** The Target data set window above is to verify the correct data set and can be switched off by enable the 'Do not ask again' option.

16. Click on 

14. Click on  to close the Shape Tool window

1. Select the '**AcquPars**' tab by clicking on it
2. Click on  to display the pulse program parameters
3. Make the following changes:


NS = 64  
DS = 8  
D1 = 2  
D9 = 0.075  
SPOFF2 = value from 7.1.5, step 14 (e.g. 645.47)  
GPZ1 = 15

1. Select the '**Spectrum**' tab by clicking on it
2. Type **zg** to start the acquisition
3. Type **ef**
4. Phase the spectrum using the manual phase adjust



NOTE: All peaks should be phased positive.

Figure 7.12.



## Plotting the reference and the TOCSY spectrum on to the same page.

7.1.6





1. Type **re 2** to display the selective TOCSY spectrum
2. Click on 
3. Type **re 1** on the command line (reference spectrum)
4. Click on 
5. Using the display tools  to adjust the spectra

Figure 7.13.



6. Type **prnt** on the command line to print the active window



NOTE: To plot the two spectra using the plot editor, follow the instructions in the manual Step-by-Step Tutorial, Basic Experiments Users Guide, Chapter 7, Homodecoupling, 7.1.6 Plotting the reference and decoupled spectra on the same page, steps 1 through 21.



# 1-D DEPT using a shaped $^{13}\text{C}$ pulse

# 8

## Introduction

8.1



Using this experiment will yield a higher Signal to noise compared with the conventional DEPT135. It is more noticeable on higher field instrument using a larger sweep width. To run this experiment the instrument has to be equipped with the hardware to do Shaped Pulses.

### Sample:

30 mg Brucine in  $\text{CDCl}_3$

## Experiment set up

8.1.1


1. Type **edc** and change the following parameters

Figure 8.1.

The screenshot shows a 'New...' dialog box with the following fields and values:

|                       |   |
|-----------------------|---|
| NAME                  | spdept  |
| EXPNO                 | 1   |
| PROCNO                | 1   |
| DIR                   | C:  |
| USER                  | pz  |
| Solvent               | $\text{CDCl}_3$   |
| Experiment Dirs.      | C:/Bruker/TOPSPIN2.1.1/exp/star/nmr/par   |
| Experiment            | C13DEPT135  |
| TITLE                 | 30 mg Brucine in $\text{CDCl}_3$<br>1D DEPT135 experiment using a 180 deg. shaped pulse in f1 |
| Receivers (1,2, ...8) | 1   |

Buttons at the bottom: OK, Cancel, More Info..., Help

2. Click on 
3. Insert the sample
4. Type **lock** and select CDCI3
5. Shim for best homogeneity
6. Tune the probe
7. Type **pulprog deptsp135** in the command line
8. Type **getprosol** in the command line

## Calculating the shaped pulse power level

8.1.2




1. In the main menu click on **'Spectrometer'** and select **'Shape Tool'** or type **stdisp** in the command line
2. In the shape tool menu bar click on  and select **'Open Shape'**

Figure 8.2.




3. Select **'Crp60comp.4'**
4. Click on 

Figure 8.3.



5. In the main menu click on **'Analysis'** and select **'Integrate Adiabatic Shape'**

Figure 8.4.



6. Make the following change:

Length of pulse [usec] = **2000**

7. Press the 'Enter' key

Figure 8.5.

**CompositeSmoothedChirp Add Const. Phase -84.0**

Integradia

2000 Length of pulse [usec]

11.25 90 deg hard pulse [usec]

Results

1.2000E08 Sweep rate on resonance [Hz/sec]

4370.19 GammaB1(max)/2pi/sqrt(Q) [Hz]

25.5832 Corresp. 90 deg square pulse [usec]

7.1360 Change of power level comp. to hard pulse [dB]

Calculator

5.0 Q (for middle of shape)

9772.05 GammaB1(max)/2pi [Hz]

update parameters



**NOTE:** The value for 'change of power lev comp. to lev of hard pulse' is calculated after executing Figure 8.5.

8. In the main menu click on '**Options**' and select '**Define Parameter Table**' by clicking on it

Figure 8.6.

|                             |        |      |
|-----------------------------|--------|------|
| Options                     | Window | Help |
| Set Path to Shape Directory |        |      |
| Define Parameter Table      |        |      |
| Select associated dataset   |        |      |
| Preferences...              |        |      |
| Remote Connection...        |        |      |
| Administration              |        |      |

Figure 8.7.

**Define Parameter Table**

Parameter Table

P 1 Length of hard pulse

PL 1 Power level of hard pulse

P 12 Length of shaped pulse

SP 2 Power level of shaped pulse

SPNAM2 Name of shaped pulse

P 13 Length of ref. shaped pulse

SP 3 Power level of ref shaped pulse

SPNAM3 Name of ref shaped pulse


FQ1LIST Name of frequency list

OK Cancel



9. Make the following changes:

- Length of hard pulse = P1
- Power level of hard pulse = PL1
- Length of shaped pulse = p12
- Power Level of shaped pulse = SP2
- Name of shaped pulse = SPNAM2

10. Click on 





11. Click on 

Figure 8.8.



12. Click on 


13. Click on  to close the Shape Tool window

### Running the experiment

8.1.3

1. Type **rga**
2. Type **zg** to start the acquisition
3. Process and Phase correct the spectrum

Figure 8.9.





# 2-D HSQC using a shaped $^{13}\text{C}$ pulse

# 9

## Introduction

9.1



Using this experiment will yield a higher Signal to noise compared with the conventional HSQCETGP. It is more noticeable on higher field instrument using a larger sweepwidth. To run this experiment the instrument has to be equipped with the hardware to do Shaped Pulses and Gradients.

### Sample:

30mg Brucine in  $\text{CDCl}_3$

## Reference spectrum

9.1.1

1. Type **edc** and change the following parameters

Figure 9.1.

The screenshot shows a 'New...' dialog box with the following fields and values:

|                       |   |
|-----------------------|---|
| NAME                  | selhsqc   |
| EXPNO                 | 1   |
| PROCNO                | 1   |
| DIR                   | C:  |
| USER                  | pz  |
| Solvent               | $\text{CDCl}_3$   |
| Experiment Dirs.      | C:/Bruker/TOPSPIN2.1.1/exp/star/nmr/par                   |
| Experiment            | PROTON  |
| TITLE                 | 1-D Proton experiment<br>30 mg Brucine in $\text{CDCl}_3$ |
| Receivers (1,2, ...8) | 1   |

Buttons at the bottom: OK, Cancel, More Info..., Help




2. Click on 
3. Insert the sample
4. Type **lock** and select  $\text{CDCl}_3$
5. Shim for best homogeneity
6. Select the '**AcquPars**' tab by clicking on it
7. Click on  to read in the Prosol parameters
8. Tune the probe
9. Type **rga**
10. Type **zg** to start the acquisition
11. Process and Phase correct the spectrum

Figure 9.2.




### Limit setting

9.1.2

1. Type **wrpa 2** on the command line
2. Type **re 2**
3. Expand the spectrum to display all peaks, leaving ca. 0.5 ppm of baseline on either side of the spectrum

Figure 9.3.






4. Click on  to set the sweep width and the O1 frequency of the displayed region

Figure 9.4.



5. Write down the SW value, rounding off to the nearest 1/10th of a ppm (e.g. **8.8**)
6. Write down the O1 value, rounding off to the nearest Hz (e.g. **1355**)
7. Click on 
8. Type **sr** on the command line and write down the exact value (e.g. **0 Hz**)

### Running the 2-D HSQC using a 180 adiabatic inversion shaped pulse in F1

9.1.3

1. Type **rpar HSQCETGPSISP all**
2. Turn the spinner off



NOTE: 2-D experiments should be run non spinning


3. Select the '**AcquPars**' tab by clicking on it
4. Make the following changes:  
SW F2 [ppm] = e.g. **8.8** (Limit setting 9.1.2, step 5)  
O1 [Hz] = e.g. **1355** (Limit setting 9.1.2, step 6)  
SOLVENT = **CDCl3**



---

All Bruker 2D inverse parameter sets use  $^{13}\text{C}$  in the F1 dimension and the sweep width and O1 are optimized to include all Carbon peaks of interest. For HSQC experiments the sw is optimized to 160 ppm.

---

5. Click on  to read in the Prosol parameters



---

The values for the pulse length and power level of the 180 deg. adiabatic inversion shaped pulse (crp60,0.5.20.1) have to be entered in to the prosol table.

---

6. Select the '**ProcPar**' tab by clicking on it
7. Make the following changes:  
SR F2 = e.g. **0** (Limit setting 9.1.2, step 8)
- 8 Select the '**Title**' tab by clicking on it
9. Change the title to: **30 mg Brucine in CDCl3, 2D HSQC using a 180 deg adiabatic inversion shaped pulse in F1**
10. Select the '**Spectrum**' tab by clicking on it

### Acquisition

9.1.4

1. Type **rga**
2. Type **zg** to start the acquisition



The standard Bruker parameter sets are optimized to run under complete automation. One of the processing parameters is an AU program for processing the data, which can be executed with the command 'xaup'. The next steps assures to use the external spectrum of Brucine for the F2 and F1 projections.

1. Type **edc2**

Figure 9.5.

A screenshot of a Windows-style dialog box titled 'edc2'. The dialog box contains a table with two columns for specifying data sets 2 and 3. The table has five rows for parameters: NAME, EXPNO, PROCNO, DIR, and USER. The first column is for data set 2 and the second for data set 3. At the bottom right of the dialog are 'OK' and 'Cancel' buttons.

| Please specify data sets 2 and 3: |         |         |
|-----------------------------------|---------|---------|
| NAME =                            | selhsqc | selhsqc |
| EXPNO =                           | 1       | 3       |
| PROCNO =                          | 1       | 1       |
| DIR =                             | c:      | c:      |
| USER =                            | pz      | pz      |



2. Enter the EXPNO and PROCNO of the 1D Proton spectrum into the first column (data set 2)
3. Click on 
4. Type **xaup**

Figure 9.6.



---

The processing AU program includes the 2D Fourier transform, baseline correction and plotting of the data.

---



# 2-D Selective HMBC

# 10

## Introduction

10.1



NOTE: To run this experiment the instrument and probe has to be equipped with the Gradient hardware.

### Sample:

50 mM Gramicidin S in DMSO

## Reference spectrum

10.1.1

1. Type **edc** and change the following parameters

Figure 10.1.

New...

Prepare for a new experiment by creating a new data set and initializing its NMR parameters according to the selected experiment type. For multi-receiver experiments several datasets are created. Please define the number of receivers in the box below.

|                       |  |
|-----------------------|--|
| NAME                  | selhmbc                                    |
| EXPNO                 | 1  |
| PROCNO                | 1  |
| DIR                   | C:   |
| USER                  | pz   |
| Solvent               | DMSO                                       |
| Experiment            | PROTON                                     |
| TITLE                 | 50mM Gramicidin S in DMSO d6<br>1-D Proton |
| Receivers (1,2, ...8) | 1  |

OK Cancel More Info... Help

2. Click on



3. Insert the sample
4. Type **lock** and select DMSO
5. Shim for best homogeneity
6. Select the '**AcquPars**' tab by clicking on it
7. Click on  to read in the Prosol parameters
8. Tune the probe
9. Type **rga**
10. Type **zg** to start the acquisition
11. Process and Phase correct the spectrum

Figure 10.2.




### Limit setting

10.1.2

1. Type **wrpa 2** on the command line
2. Type **re 2**
3. Expand the spectrum to include all peaks (e.g. 0 ppm to 10 ppm)

Figure 10.3.






4. Click on  to set the sweep width and the O1 frequency of the displayed region

Figure 10.4.



5. Click on 
6. Type **sw** on the command line and write down the value of SW, rounding off to the nearest 1/10th of a ppm (e.g. **10 ppm**)
7. Type **o1p** on the command line and write down the value of O1p, rounding off to the nearest 1/10th of a ppm (e.g. **5 ppm**)
8. Type **sr** on the command line and write down the exact value (e.g. **0 Hz**)

### Running a 2-D HMBC experiment

10.1.3

1. Type **rpar HMBCGP all**

2. Turn the spinner off



---

NOTE: 2-D experiments should be run non spinning

---

3. Select the '**AcquPars**' tab by clicking on it

4. Make the following changes:

F2 SW [ppm] = value from step 6, Limit setting 10.1.2 (e.g. **10**)

O1P [Hz] = value from step 7, Limit setting 10.1.2 (e.g. **5**)


Solvent = DMSO



---

All Bruker 2D inverse parameter sets use <sup>13</sup>C in the F1 dimension and the sweep width and O1 are optimized to include all Carbon peaks of interest. For HMBC experiments the sw is optimized to 220 ppm.

---

5. Click on  to read in the Prosol parameters

6. Select the '**ProcPar**' tab by clicking on it

7. Make the following changes:

SR [F2] = value from step 8, Limit setting 10.1.2 (e.g. **0**)

8 Select the '**Title**' tab by clicking on it

9. Change the title to: **50 mM Gamicidin S in DMSO-d6, 2-D HMBC**

10. Select the '**Spectrum**' tab by clicking on it

---

### Acquisition

10.1.4

1. Type **rga**

2. Type **zg** to start the acquisition



The standard Bruker parameter sets are optimized to run under complete automation. One of the processing parameters is an AU program for processing the data, which can be executed with the command 'xaup'. The next steps assures to use the external spectrum of Gramicidin for the F2 projection.

1. Type **edc2**

Figure 10.5.

| Please specify data sets 2 and 3: |         |         |
|-----------------------------------|---------|---------|
| NAME =                            | selhmbc | selhmbc |
| EXPNO =                           | 1       | 2       |
| PROCNO =                          | 1       | 3       |
| DIR =                             | C:      | C:      |
| USER =                            | pz      | pz      |

2. Enter the EXPNO and PROCNO of the 1D Proton spectrum into the first column (data set 2)


3. Click on



The processing AU program includes the 2D Fourier transform, baseline correction and plotting of the data. The HMBC experiment uses magnitude mode for processing and shows only positive peaks.

4. Type **xaup**

Figure 10.6.



### Optimizing the parameters on the carbonyl region

10.1.6

1. Type **wrpa 3** on the command line
2. Type **re 3**
3. Expand the carbonyl region including all cross peaks (e.g. 168 ppm to 178 ppm)

Figure 10.7.






Figure 10.8.




4. Write down the expanded F1 sweep width in ppm (e.g. **10 ppm**)
5. Write down the center frequency (O2) of the expanded F1 sweep width in ppm (e.g. **172 ppm**)
6. Select the '**AcqPars**' tab by clicking on it
7. Click on  to display the pulse program parameters
8. Write down the value for P3 [us] (e.g. **9.5 us**)
9. Write down the value for PL2 [dB] (e.g. **0 dB**)
10. Select the '**Title**' tab by clicking on it
11. Change the title to: **50 mM Gamicidin S in DMSO-d6, selective 2-D HMBC**
12. Select the '**Spectrum**' tab by clicking on it

### Set up the selective pulse

10.1.7


1. Type **pulprog shmbcgpnd** in the command line
2. In the main menu click on '**Spectrometer**' and select '**Shape Tool**' or type **stdisp** in the command line

Figure 10.9.



3. In the main menu click on 'Shapes', select 'Classical' and select 'Sinc' by clicking on it

Figure 10.10.




4. Make the following changes:


Change size of shape = 256

Number of cycles = 3



Figure 10.11.



5. Click on 
6. Click on **'Save Shape'**
7. Make the following changes:

File Name = **Sinc3.256**

Figure 10.12.





8. Click on 
9. In the main menu click on **'Analysis'**, select **'Calculate Bandwidth for Excitation'**

Figure 10.13.



10. Make the following changes:


DeltaOmega [Hz] = **750** (e.g. SW 10 ppm from step 4 in 10.1.6)

11. Press the 'Enter' key



NOTE: The value of Delta T [usec] is being calculated. (e.g. 7450.7 usec)


Figure 10.14.



12. Write down the Delta T value [usec] (e.g. **7450.7 usec**)

13. Click on


Figure 10.15.



14. Click on

15. In the main menu click on '**Analysis**', select '**Integrate Shape**'

Figure 10.16.



16. Make the following change:

Length of pulse [usec] = value from step 12 in 10.1.7 (e.g. **7450.7**)

17. Press the 'Enter' key


Total rotation [degree] = **90**

18. Press the 'Enter' key


90 deg. hard pulse [usec] = p3 from step 8 in 10.1.6 e.g. **9.5**)

19. Press the 'Enter' key

Figure 10.17.



20. Write down the change of power level [dB] value (e.g. **42.85**)

21. Click on  to close the Shape Tool window

## Setting up the acquisition parameters

10.1.8


1. Select the 'AcquPars' tab by clicking on it

2. Make the following changes:

NS = 32

F1 SW [ppm] = value from step 4 in 10.1.6 (e.g. 10)

O2P [ppm] = value from step 5 in 10.1.6 (e.g. 172)

3. Click on  to display the pulse program parameters

4. Make the following changes:

P13 [us] = value from step 12 in 10.1.7 (e.g. 7450.7)


SP14 [dB] = (value from step 19 in 10.1.7) + (PL2) (e.g. 42.85)

## Running the experiment

10.1.9


1. Select the 'Spectrum' tab by clicking on it
2. Type **zg** to start the acquisition
3. Type **xfb** to process the 2-D data
4. Expand the 2-D spectrum

Figure 10.18.



5. Compare the result of the selective HMBC against the regular HMBC in 10.1.3

Figure 10.19.



---

NOTE: The cross peaks in the selective HMBC show nice separation do to the increased resolution in F1, compared to the regular HMBC. The projections are external high resolution spectra.

---





## Notes: